

# C-Series Rotary Valve

## FEATURES

- > 2-year warranty against defects in materials and workmanship
- > Long-life
- > Exceptional Value
- > Available in a variety of inert materials
- > No solenoid valve limitations (leaking, clogging, low CVs)
- > No ceramic valve limitations (seizing)

## CONFIGURED OPTIONS

- > Materials
  - Body PEEK  
PCTFE
  - Rotary Plug UHMWPE  
PTFE
- > Valve Ports
  - "Y" 120° 3-port
  - 90° 4-port
  - Distribution 4-port
  - T 3-port
  - T 4-port
  - Distribution 6-way

## OEM OPTIONS






- > Custom configurations possible

## BASE MODEL

Compatible with C-Series Pumps and Rotary Valve Controllers



**SPECIFICATIONS**

Description	Materials (Body/Plug)	Catalog Number	Flow Path (Blue)
"Y" 120° 3-port	PEEK/UHMWPE	10020	
	PEEK/PTFE	10021	
	PCTFE/UHMWPE	10022	
	PCTFE/PTFE	10023	
90° 4-port	PEEK/UHMWPE	10040	
	PEEK/PTFE	10041	
	PCTFE/UHMWPE	10042	
	PCTFE/PTFE	10043	
Distribution 4-Port	PCTFE/UHMWPE	10062	
	PCTFE/PTFE	10063	
T 3-port	PEEK/UHMWPE	10115	
T 4-port	PEEK/UHMWPE	10125	
Distribution 6-way	PCTFE/UHMWPE	10106	

ROTARY VALVES FOR SYRINGE PUMPS

[tricontinent.com](http://tricontinent.com)



TriContinent Scientific, Inc.  
12555 Loma Rica Drive  
Grass Valley, CA 95945  
USA  
Tel: 800-937-4738  
Fax: 530-273-2586  
[liquidhandling.tcs@gardnerdenver.com](mailto:liquidhandling.tcs@gardnerdenver.com)

Gardner Denver Thomas GmbH  
Livry-Gargan-Str. 10  
82256 Fürstentfeldbruck  
Germany  
Tel: +49 8141 2280 0  
Fax: +49 8141 8892136  
[thomas.de@gardnerdenver.com](mailto:thomas.de@gardnerdenver.com)

Gardner Denver Hong Kong, Limited  
Unit 1317-1318 Delta House, 3 On Yiu Street  
Siu Lek Yuen, Shatin, New Territories  
Hong Kong  
Tel: +852 26903502  
Fax: +852 27924598  
[thomas.hk@gardnerdenver.com](mailto:thomas.hk@gardnerdenver.com)



The information presented in this material is based on technical data and test results of nominal units. It is believed to be accurate and reliable and is offered as an aid to help in the selection of TriContinent products. It is the responsibility of the user to determine the suitability of the product for the intended use and the user assumes all risk and liability in connection there with. TriContinent does not warrant, guarantee or assume any obligation or liability in connection with this information.

Photos of products pictured in this catalog do not necessarily represent a specific model number. To obtain further information for custom options, contact your local TriContinent office.

Printed in USA Form No. MKT90034 C © Gardner Denver TriContinent. All rights reserved.